

Genliste Archer Myeloid Panel (TruSight® Myeloid Sequencing Panel)
Mutationen

ABL1	ASXL1	ATRX	BCOR	BCORL1	BRAF
Exon 4-10		Exon 8-11, 17-32	Exon 2-15		Exon 3, 10-13 und 15
CALR	CBL	CBLB	CBLC	CDKN2A	CEBPA
Exon 8 und 9	Exon 2-5, 7-9 und 16	Exon 3, 9 und 10	Exon 9 und 10		
CSF3R	CUX1	DDX41	DNMT3A	ETNK1	ETV6
Exon 10 und 14-18				Exon 3	
EZH2	FBXW7	FLT3	GATA1	GATA2	GNAS
Exon 2-20		Exon 8-17, 19-21	Exon 2	Exon 2-6	Exon 8-11
HRAS	IDH1	IDH2	IKZF1	JAK2	JAK3
Exon 2-4	Exon 3 und 4	Exon 4 und 6		Exon 12-16 und 19-25	Exon 3,11,13,15,18,19
KDM6A	KIT	KRAS	MLL/KMT2A	MPL	MYD88
	Exon 1,2,5,8-15,17,18	Exon 2-4		Exon 10 und 12	Exon 3-5
NF1	NOTCH1	NPM1	NRAS	PDGFRA	PHF6
	Exon 26-28, 34	Exon 12	Exon 2-5	Exon 12, 14, 15 und 18	Exon 2-10
PPM1D	PTEN	PTPN11	RAD21	RUNX1	SETBP1
Exon 6		Exon 3, 4, 7, 8, 12, 13	Exon 2-14		Exon 4
SF3B1	SH2B3	SMC1A	SMC3	SRSF2	STAG2
Exon 13-21	Exon 2-8		Exon 10,13,19,23,25,28	Exon 1 und 2	Exon 2-33
TET2	TP53	U2AF1	U2AF2	WT1	ZRSR2
Exon 4-11		Exon 2, 6 und 7			